

Obtaining 1D-INADEQUATE Spectra

A 1D-INADEQUATE spectrum will show carbon-carbon coupling constants. Note that this experiment filters the “main” ^{13}C peaks away, so only the sidebands due to adjacent ^{13}C atoms appear. The procedure is therefore *slow*. A concentrated sample (which will provide a 4-scan carbon spectrum) is essential; count on acquiring for a full day.

1. Acquire a normal ^{13}C spectrum:

```
rpar C13CPD all
atma
getprosol
ns=4, ds=0
```

2. Save this to a new experiment number (iexpno) and convert it to an INADEQUATE parameter set.

```
pulprog = inadgpspid
d4      = 7.6 ms            $\frac{1}{4}J_{CC}$ 
cnst3   = 33              $J_{CC}$  — used to calculate d4?
ns      = something huge  999,999?
ds      = 0? 1?
```

Type gpz and set

```
gpz1    = 40             gradient power (% of full)
gpnam1  = SMSQ10.100
gpz2    = 40
gpnam2  = SMSQ10.100
```

3. zg as usual.

Let this run for a bit, then check it:

```
tr efp    (transfer data to working memory, then efp)
```

Most likely there will be nothing to see.

Let it run a few more hours and check again. Let run overnight and check again.